

Type of test	Specimen sample	Test	Time (min)	Sensitivity (%)	Specificity (%)	Cross-reactivity
NAAT	Nasopharyngeal swab, sputum, bronchoalveolar lavage fluid	rRT-PCR	240	71–98	95	No
		CRISPR	40	97	100	No
		LAMP ^b	30–60	75–90	93–99	No
Antigen detection	Nasopharyngeal swab	-	15–30	62–92	100	Yes
Immune assay	Blood	-	15–30	92–100	93–100	Yes